# **ECG Anomaly Classification**

```
In [2]: # Import packages
        import glob
        import random
        from collections import OrderedDict
        from biosppy.signals import ecg
        import pandas as pd
        import numpy as np
        import matplotlib.pyplot as plt
        import seaborn as sns
        import scipy.interpolate as interp
        from sklearn.preprocessing import MinMaxScaler, StandardScaler
        from joblib import Parallel, delayed
        import tensorflow as tf
        from keras.models import Model, Sequential
        from keras.layers import Input, Dense, LSTM, Dropout
        from keras.preprocessing import sequence
        from keras.backend.tensorflow_backend import set_session
        from sklearn.metrics import classification_report, confusion_matrix
        from sklearn.utils import shuffle
        seed = 9441
        random.seed(seed)
        np.random.seed(seed)
        config = tf.ConfigProto()
        config.gpu_options.allow_growth = True
        config.gpu options.visible device list = "0"
        set session(tf.Session(config=config))
```

Using TensorFlow backend.

## Working on the preprocessed dataset

There are 4 types of anomalies: Control, High P Amplitude, SA and ST.

All data under 1000 Hz sampling frequency.

```
In [3]: con_dir = 'cleaned_data/control'
highp_dir = 'cleaned_data/high p wave'
sa_dir = 'cleaned_data/sinoatrial arrest'
st_dir = 'PROCESSED ECG database/processedST'
sampling_freq = 1000
In [4]: def chunks(1, n):
```

```
In [4]: def chunks(1, n):
    """Return successive n-sized chunks from l."""
    chunked_list = []
    for i in range(0, len(1), n):
        sublist = l[i:i + n]
        if len(sublist) == n:
            chunked_list.append(l[i:i + n])
    return chunked_list
```

```
In [5]: def get_len_sublists(alist, level=1):
    if level == 1:
        return [len(sublist) for sublist in alist]
    elif level == 2:
        lens = []
        for inlist in alist:
            lens.extend(get_len_sublists(inlist))
        return lens
```


```
In [6]: def select_random(alist):
    aux = list(alist)
    np.random.shuffle(aux)
    return aux[0]
```

Original signals whose lengths are greater than 10k data points (10 seconds of data recorded) are split into multiple samples of the same length of 10k if possible.

```
In [7]: def read signals(adir, all files=False):
            # If a signal sequence is too long, split it into sublists
             standard_size = 10000 # 10 seconds
            if all_files:
                files = sorted(glob.glob(adir + '/*.csv'))
            else:
                files = sorted(glob.glob(adir + '/ *.csv'))
            signals = [pd.read_csv(fi, header=None)[0].values for fi in files]
            standard_signals = [signal for signal in signals if len(signal) < standard_size]</pre>
            oversized_signals = [signal for signal in signals if len(signal) >= standard_size]
            for signal in oversized_signals:
                 standard_signals.extend(chunks(signal, standard_size))
            return standard signals
        con_signals = read_signals(con_dir, all_files=True)
        highp_signals = read_signals(highp_dir, all_files=True)
        sa_signals = read_signals(sa_dir, all_files=True)
        st_signals = read_signals(st_dir)
```

```
In [8]: anomalies = ['Control', 'High P', 'SA', 'ST']
    y_pos = np.arange(len(anomalies))
    lengths = [len(con_signals), len(highp_signals), len(sa_signals), len(st_signals)]


plt.bar(y_pos, lengths, align='center')
    plt.xticks(y_pos, anomalies)
    plt.ylabel('No of Signals')
    plt.title('Total signals')
    plt.show()
```




```
In [9]: plt.rcParams['figure.figsize'] = [15, 7]
```


# Plot a Control ECG Singal







# Plot a High P Amplitude ECG Singal







Plot a SA ECG Singal





#### Plot a ST ECG Singal






In [14]: plt.rcParams['figure.figsize'] = [6, 4]

### Scale signal amplitudes


#### **Extract R-R intervals**


```
In [16]: def cal r peaks(signal, sampling rate=sampling freq):
             rpeaks, = ecg.hamilton segmenter(signal=signal, sampling rate=sampling rate)
             rpeaks, = ecg.correct_rpeaks(signal=signal, rpeaks=rpeaks, sampling_rate=sampling_rate, tol=0.05)
             templates, rpeaks = ecg.extract_heartbeats(signal=signal, rpeaks=rpeaks, sampling rate=sampling_ra
         te, before=0.2, after=0.4)
             return rpeaks
         def get_r_peaks(signals):
             return Parallel(n_jobs=6)(delayed(cal_r_peaks)(signal) for signal in signals)
In [17]:
         con r peaks = get r peaks(con signals scaled)
         highp r peaks = get r peaks(highp signals scaled)
         sa_r peaks = get r_peaks(sa signals scaled)
         st_r_peaks = get_r_peaks(st_signals_scaled)
In [18]: def split_list(alist, indices):
             splitted list = []
             for start, end in zip(indices, indices[1:]):
                 sublist = alist[start:end+1]
                 splitted_list.append(sublist)
             return splitted_list
         def extract rr intervals(signals, r peaks):
             rr itvs = []
             for idx, signal in enumerate(signals):
                 if r_peaks[idx] is None:
                     continue
                 rr_itvs.append(split_list(signal, r_peaks[idx].tolist()))
             return rr itvs
         con rr itvs = extract rr intervals(con signals scaled, con r peaks)
         highp_rr_itvs = extract_rr_intervals(highp_signals_scaled, highp_r_peaks)
         sa_rr_itvs = extract_rr_intervals(sa_signals_scaled, sa_r_peaks)
         st_rr_itvs = extract_rr_intervals(st_signals_scaled, st_r_peaks)
In [20]: plt.rcParams['figure.figsize'] = [6, 4]
         def plot_rr_intervals(rr_intervals, anomaly):
             for rr in rr_intervals[1:-1]:
                 plt.plot(rr)
             plt.title(anomaly)
             plt.show()
```

In [21]: plot\_rr\_intervals(select\_random(con\_rr\_itvs), 'Con')
 plot\_rr\_intervals(select\_random(highp\_rr\_itvs), 'HighP')
 plot\_rr\_intervals(select\_random(sa\_rr\_itvs), 'SA')
 plot\_rr\_intervals(select\_random(st\_rr\_itvs), 'ST')



### R-R intervals are normalized to a specific length





```
In [24]:
         rr normalized size = 1000
         max_rr_intervals = 18
         def linear_interpolation(rr_interval):
             rr_interp = interp.interp1d(np.arange(rr_interval.size), rr_interval)
             return rr_interp(np.linspace(0,rr_interval.size-1, rr_normalized_size))
         def pre pad sequence(sequence):
             max_length = max_rr_intervals
             seq_length = sequence.shape[0]
             if seq_length == max_length:
                 return sequence
             elif seq_length < max_length:</pre>
                 dim = max length - seq length
                 zeros seq = np.zeros((max length - seq length, rr_normalized_size))
                 return np.concatenate([zeros_seq, sequence])
             else:
                 return sequence[seq_length - max_length:]
         def normalize_rr_intervals(signal_rr_intervals):
             signal rr itvs normed = []
             for signal in signal_rr_intervals:
                 rr_itvs_normed = []
                 if len(signal) == 0:
                     continue
                 for rr in signal:
                     rr_itvs_normed.append(linear_interpolation(rr))
                 rr itvs normed padded = pre pad sequence(np.array(rr itvs normed))
                 signal rr itvs normed.append(np.array(rr itvs normed padded))
             return np.array(signal_rr_itvs_normed)
         con_rr_itvs_normed = normalize_rr_intervals(con_rr_itvs)
         highp_rr_itvs_normed = normalize_rr_intervals(highp_rr_itvs)
         sa_rr_itvs normed = normalize_rr_intervals(sa_rr_itvs)
         st_rr_itvs_normed = normalize_rr_intervals(st_rr_itvs)
```

## R-R Interval Dimension Reduction with Autoencoder

```
In [27]:
         #Split train, evaluation, and test sets
         eval_split_pos = int(len(y) * .7)
         test_split_pos = int(len(y) * .8)
         X_train, y_train = X[:eval_split_pos], y[:eval_split_pos]
         print('X_train.shape', X_train.shape)
         X_eval, y_eval = X[eval_split_pos:test_split_pos], y[eval_split_pos:test_split_pos]
         print('X_eval.shape', X_eval.shape)
         X_test, y_test = X[test_split_pos:], y[test_split_pos:]
         print('X_test.shape', X_test.shape)
         # Reshape rr sequences into unordered rr intervals and remove zero padded rr rows
         X_train_AE = X_train.reshape(-1, rr_normalized_size)
         print('X_train.reshape', X_train_AE.shape)
         all zeros rows train = (X train AE==0).all(1)
         X train AE = X train AE[~all zeros rows train]
         print('X_train_AE.shape', X_train_AE.shape)
         X eval AE = X eval.reshape(-1, rr normalized size)
         print('X_eval.reshape', X_eval_AE.shape)
         all zeros rows eval = (X eval AE==0).all(1)
         X eval AE = X eval AE[~all zeros rows eval]
         print('X eval AE.shape', X eval AE.shape)
         X_test_AE = X_test.reshape(-1, rr_normalized_size)
         print('X_test.reshape', X_test_AE.shape)
         all_zeros_rows_test = (X_test_AE==0).all(1)
         X_test_AE = X_test_AE[~all_zeros_rows_test]
         print('X_test_AE.shape', X_test_AE.shape)
         ('X_train.shape', (4506, 18, 1000))
         ('X_eval.shape', (644, 18, 1000))
         ('X_test.shape', (1288, 18, 1000))
         ('X_train.reshape', (81108, 1000))
         ('X_train_AE.shape', (30079, 1000))
         ('X_eval.reshape', (11592, 1000))
         ('X_eval_AE.shape', (4267, 1000))
         ('X_test.reshape', (23184, 1000))
('X_test_AE.shape', (8694, 1000))
```

```
Train on 30079 samples, validate on 4267 samples
Epoch 1/25
Fnoch 2/25
Epoch 3/25
Epoch 4/25
Epoch 5/25
Fnoch 6/25
Epoch 7/25
Epoch 8/25
Epoch 9/25
Epoch 10/25
Epoch 11/25
Epoch 12/25
Epoch 13/25
Epoch 14/25
Epoch 15/25
Fnoch 16/25
Epoch 17/25
Epoch 18/25
Epoch 19/25
Fnoch 20/25
Epoch 21/25
Epoch 22/25
Epoch 23/25
Epoch 24/25
Epoch 25/25
```

http://213.246.38.101:9441/ 10/14

Out[28]: <keras.callbacks.History at 0x7fa47c431890>

#Indices to reinsert zero padded rows

In [29]:

```
from operator import itemgetter
         from itertools import groupby
         def __get_continuous_ranges(indices):
              ranges = []
              for k, g in groupby(enumerate(indices), lambda (i,x):i-x):
                  group = map(itemgetter(1), g)
                  ranges.append((group[0], group[-1]))
              return ranges
         zeros indices train = np.nonzero(all zeros rows train)[0]
         zeros indices eval = np.nonzero(all zeros rows eval)[0]
         zeros indices test = np.nonzero(all zeros rows test)[0]
         zero_ranges_train = __get_continuous_ranges(zeros_indices_train)
         zero_ranges_eval = __get_continuous_ranges(zeros_indices_eval)
zero_ranges_test = __get_continuous_ranges(zeros_indices_test)
In [30]: #Encode rr intervals
         X_train_AE_encoded = encoder.predict(X_train_AE)
         print('X train AE encoded.shape', X train AE encoded.shape)
         X_eval_AE_encoded = encoder.predict(X_eval_AE)
         print('X_eval_AE_encoded.shape', X_eval_AE_encoded.shape)
         X test AE encoded = encoder.predict(X test AE)
         print('X_test_AE_encoded.shape', X_test_AE_encoded.shape)
         #Reinsert zero padded rows
         for arange in zero ranges train:
              X train AE encoded = np.insert(X train AE encoded, arange[0], np.zeros((arange[1] - arange[0] + 1,
          encoding dim)), 0)
         print('X train AE encoded.shape', X train AE encoded.shape)
         for arange in zero ranges eval:
              X eval AE encoded = np.insert(X eval AE encoded, arange[0], np.zeros((arange[1] - arange[0] + 1, e
         ncoding_dim)), 0)
         print('X_eval_AE_encoded.shape', X_eval_AE_encoded.shape)
         for arange in zero ranges test:
              X test AE encoded = np.insert(X test AE encoded, arange[0], np.zeros((arange[1] - arange[0] + 1, e
         ncoding dim)), 0)
         print('X_test_AE_encoded.shape', X_test_AE_encoded.shape)
         # Reshape to rr sequences
         X_train_dim_reduced = X_train_AE_encoded.reshape((-1, max_rr_intervals, encoding_dim))
         print('X_train_dim_reduced.shape', X_train_dim_reduced.shape)
         X eval dim reduced = X eval AE encoded.reshape(-1, max rr intervals, encoding dim)
         print('X eval dim reduced.shape', X eval dim reduced.shape)
         X_test_dim_reduced = X_test_AE_encoded.reshape(-1,max_rr_intervals, encoding_dim)
         print('X test dim reduced.shape', X test dim reduced.shape)
         ('X_train_AE_encoded.shape', (30079, 500))
         ('X_eval_AE_encoded.shape', (4267, 500))
         ('X_test_AE_encoded.shape', (8694, 500))
         ('X_train_AE_encoded.shape', (81108, 500))
         ('X_eval_AE_encoded.shape', (11592, 500))
         ('X_test_AE_encoded.shape', (23184, 500))
         ('X_train_dim_reduced.shape', (4506, 18, 500))
('X_eval_dim_reduced.shape', (644, 18, 500))
         ('X test_dim_reduced.shape', (1288, 18, 500))
```

### **Anomaly Classification with LSTM**

```
Layer (type)
            Output Shape
                        Param #
-----
1stm 1 (LSTM)
            (None, 64)
                        144640
dropout 1 (Dropout)
            (None, 64)
                        a
dense 5 (Dense)
                        260
            (None, 4)
Total params: 144,900
Trainable params: 144,900
Non-trainable params: 0
None
Train on 4506 samples, validate on 644 samples
Epoch 1/20
0.5085 - val_acc: 0.8137
Epoch 2/20
0.4174 - val_acc: 0.8478
Epoch 3/20
0.4607 - val acc: 0.8106
Epoch 4/20
4506/4506 [==================== ] - 3s 613us/step - loss: 0.3992 - acc: 0.8844 - val loss:
0.3052 - val acc: 0.8944
Epoch 5/20
0.2312 - val acc: 0.9193
Epoch 6/20
4506/4506 [=============== ] - 3s 613us/step - loss: 0.2350 - acc: 0.9292 - val loss:
0.1820 - val acc: 0.9270
Epoch 7/20
0.1737 - val_acc: 0.9224
Epoch 8/20
4506/4506 [============== ] - 3s 615us/step - loss: 0.1606 - acc: 0.9527 - val loss:
0.1042 - val_acc: 0.9627
Epoch 9/20
0.0992 - val_acc: 0.9674
Epoch 10/20
0.1681 - val acc: 0.9379
Epoch 11/20
0.0906 - val acc: 0.9627
Epoch 12/20
4506/4506 [============== ] - 3s 618us/step - loss: 0.1214 - acc: 0.9627 - val loss:
0.0607 - val_acc: 0.9860
Epoch 13/20
0.0572 - val_acc: 0.9783
Epoch 14/20
0.0453 - val_acc: 0.9829
Epoch 15/20
0.0392 - val acc: 0.9860
Epoch 16/20
0.0732 - val_acc: 0.9689
Epoch 17/20
0.1497 - val acc: 0.9317
Epoch 18/20
0.1025 - val_acc: 0.9705
Epoch 19/20
```

```
0.0578 - val_acc: 0.9876
        Epoch 20/20
        0.0350 - val_acc: 0.9938
Out[31]: <keras.callbacks.History at 0x7fa4a014f610>
In [32]: y_preds = np.argmax(model.predict(X_test_dim_reduced), axis=1)
        y_trues = np.argmax(y_test, axis=1)
        print(classification_report(y_trues, y_preds, target_names=anomalies))
        print('Confusion Matrix:\n')
        print(confusion_matrix(y_trues, y_preds, labels=[0, 1, 2, 3]))
                   precision
                             recall f1-score support
                       0.99
           Control
                               0.99
                                        0.99
                                                  77
            High P
                       1.00
                               0.94
                                       0.97
                                                 660
                               0.99
               SA
                       0.93
                                       0.96
                                                 537
                       1.00
               ST
                               0.86
                                       0.92
                                                 14
                       0.97
                               0.97
                                       0.97
                                                1288
        avg / total
```

#### Confusion Matrix:

```
[[ 76  0  1  0]
 [ 0 622  38  0]
 [ 0  3 534  0]
 [ 1  0  1 12]]
```